

XpressCheck His Tag rapid test kit



(REF# XC-HIS, Store at Room Temperature, keep dry)

XpressCheck His Tag is a competitive lateral flow assay for rapid detection of His-tagged proteins. It provides results in 3 minutes without instrumentation and works directly with unpurified samples, including cell lysates and culture media. The strip supports both qualitative and semi-quantitative analysis, making it ideal for protein expression screening and optimization. Its detection limit is approximately 0.1 μM (about 5 $\mu\text{g}/\text{mL}$ for a 50 kDa protein).

Kit contents

- Test strips (in desiccant tube)
- Xpress buffer
- Instruction card

Result Interpretation

Result	Negative	Positive
Example	 <p>Both C and T lines are clearly visible and have similar intensity.</p>	 <p>T line is significantly weaker than C, or only C line is visible.</p>

Troubleshooting

Problem	Suggestion
Missing C line	The result is invalid, likely due to insufficient sample volume or high sample viscosity. Increase the sample volume or dilute the sample with water or buffer and repeat the test.
Does not detect old purified His-tagged proteins.	This kit is intended for screening unpurified samples or analyzing freshly purified proteins. Prolonged storage may cause purified His-tagged proteins to degrade or aggregate, potentially masking the His tag and reducing detection.

See reverse side for Application Instructions

Application Example 1: Rapidly assessing secreted recombinant protein expression in medium

Transfer 100 μ L medium to a well and place a test strip in the well.
(Dilute the sample with water or buffer if necessary)



Read result after 3 minutes

Application Example 2: Rapidly assessing recombinant protein expression in E. coli lysate

E.Coli pellet from
1 mL culture



Lyse cells with 100 μ L
Xpress buffer



Transfer to a well and place a
test strip in the well



Read result after 3 minutes

Application Example 3: Evaluation of recombinant protein expression in soluble VS inclusion body fractions

E.Coli pellet from
1 mL culture



Lyse cells following
user's protocol



Insoluble
fraction



Lyse cells with 100 μ L
Xpress buffer



Transfer to a well and place a
test strip in the well



Read result after 3 minutes

Transfer to a well and place a
test strip in the well



Read result after 3 minutes


Application Example 4: Semi-quantitative measurement of His-tagged protein concentration

Prepare serial
diluted samples


Serial
dilution



Place a test strip
in each well



Read result
after 3 minutes



The highest dilution that still produced a positive
result corresponded to approximately 0.1 μ M